

Datasheet

FAM127A mouse monoclonal antibody (hybridoma)

Catalog Number: H00008933-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant FAM127A.

Immunogen: FAM127A (NP_003919.1, 1 a.a. ~ 209 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MGGGRGLLGRET LGPGGGCSGEGPLCYWPPPGSPP
APSLRASLPLEPPRCPLRSCSLPRSACLC SRNSAPGS
CCRPWASLWSEPPSPSSQPAPPMYIWLSCAPAAS
WAPVTHWTDHPLPLPSLLPTRLPDDYIILPTDLRCH
SHRHPSHPTDRLLLLVIWTHLGGIWAGHSPWTVIQT
A
GRPPRDLSPSARPISSPPPETSCVLA

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 8933

Gene Symbol: FAM127A

Gene Alias: CXX1, MAR8C, MART8C, MGC117411, Mar8, Mart8