

Datasheet

FEM1B monoclonal antibody (M09), clone 4B12

Catalog Number: H00010116-M09

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant FEM1B.

Clone Name: 4B12

Immunogen: FEM1B (NP_056137.1, 401 a.a. ~ 490
a.a) partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

TVKAPDIECVLRCSVLEIEQSMNRVKNISDADVHNAMD
NYECNLYTFLYLVCISTKTQCSEEDQCKINKQIYNLIHL
DPRTRREGFTLLHL

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10116

Gene Symbol: FEM1B

Gene Alias: DKFZp451E0710, FIAA