

Datasheet

ZNF263 monoclonal antibody (M05), clone 2G6

Catalog Number: H00010127-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant ZNF263.

Clone Name: 2G6

Immunogen: ZNF263 (NP_005732, 201 a.a. ~ 298 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

PEGNMEDKEMTGPQLPESLEDVAMYISQEEWGHQDP
SKRALSRDVTQESYENVDSLESHIPSQEVPGTQVGQG
GKLWDPSVQSCKEGLSPRGPAPGEE

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10127

Gene Symbol: ZNF263

Gene Alias: FPM315, ZKSCAN12