

Datasheet

C4orf6 monoclonal antibody (M01), clone 3F12

Catalog Number: H00010141-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant C4orf6.

Clone Name: 3F12

Immunogen: C4orf6 (NP_005741, 1 a.a. ~ 93 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MDTQKQIHKTHNSKNQFFTIFFFLSVEFGKEGTRKNFY
LLLSIGHYGRKSRRADLGTADTADKTEPECFAASWTF
DPNPSVTVSGAHSTAVHQ

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10141

Gene Symbol: C4orf6

Gene Alias: aC1

Gene Summary: This gene is expressed in
neuroblastoma; however, the function of this gene is not
yet determined. [provided by RefSeq]