

Datasheet

CEBPZ MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00010153-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human CEBPZ protein.

Immunogen: CEBPZ (NP_005751.2, 1 a.a. ~ 1054 a.a) full-length human protein.

Sequence:

MAAVKEPLEFHAKRPWRPEEAVEDPDEEDEDNTSEA
ENGFSLEEVLRLGGTKQDYLMATLDENEEVIDGGKK
GAIDDLQQGELEAFIQNLNLAKYTKASLVEEDEPAEKE
NSSKKEVKIPKINNKNNTAESQRTSVNKVKKNRPEPHS
DENGSTTPKVKKDKQNIFEFFERQTLRLPGGKWYDL
EYSNEYSKPKPQDVVSKYKTLAQKLYQHEINLFKSKT
NSQKGASSTWMKAIVSSGTLGDRMAAMILLIQDDAVH
TLQFVETLVNLVKKKGSKQQCLMALDTFKELLITDLLPD
NRKLRIFSQRPFDKLEQLSSGNKDSRDRRLILWYFEH
QLKHLVAEFVQVLETLSHDTLVTTKTRALTVAHELLCN
KPEEEKALLVQVVNKLGDQPQNRIATKASHLLETLLCKH
PNMKGVVSGEVERLLFRSNISSKAQYYAICFLNQMAL
HEESELANKLITVYFCFFRTCVKKKDVESKMLSALLTG
VNRAYPYSQTGDDKVREQIDTLFKVLHIVNFNTSVQAL
MLLFQVMNSQQTISDRYYTALYRKMLDPGLMTCSKQA
MFLNLVYKSLKADIVLRRVKAFVKRLLQVTCQQMPPFI
CGALYLVSEILKAKPGLRSQLDHHPESDDEENFIDAND
DEDMEKFTDADKETEIVKKLETEETVPETDVETKKPEV
ASWVHFDNLKGGKQLNKYDPFSRNPLFCGAENTSLW
ELKKLSVHFHPSVALFAKTILQGNIIQYSGDPLQDFTL
MRFLDRFVYRNPKPHKGKENTDSVVMQPKRKHFIKDI
RHLPVNSKEFLAKEESQIPVDEVFFHRYKKVAVKEKQ
KRDADDEESIEDVDDEEFELIDTFEDDNCFSGKDDM
DFAGNVKKRTKGAKDNTLDEDESGSDDELGNLDDDE
VSLGSMDDDEFAEVEDDGGTFMDVLDDESESVPELE
VHSKVSTKSKRKGTDDDFAGSFQGPRKKRNLND
SSLFVSAEEFGHLLDENMGSKFDNIGMNAMANKDNAS
LKQLRWEAERDDWLHNRDAKSIKKKKHFKKKRIKTTQ
KTKKQRK

Host: Mouse

Reactivity: Human, Mouse, Rat

Applications: IF, WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10153

Gene Symbol: CEBPZ

Gene Alias: CBF2, HSP-CBF, NOC1