

Datasheet

PLXNC1 monoclonal antibody (M06), clone 1A12

Catalog Number: H00010154-M06

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant PLXNC1.

Clone Name: 1A12

Immunogen: PLXNC1 (NP_005752, 250 a.a. ~ 348 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

PYNYTSGAATGWPSMARIAQSTEVLFQQQASLDCGH
GHPDGRRLLSSSLVEALDVWAGVFSAAAGEGQERR
SPTTTALCLFRMSEIQARAKRVSWDFK

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10154

Gene Symbol: PLXNC1

Gene Alias: CD232, PLXN-C1, VESPR