

Datasheet

TRIM28 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00010155-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TRIM28 protein.

Immunogen: TRIM28 (NP_005753.1, 1 a.a. ~ 835 a.a) full-length human protein.

Sequence:

MAASAAAASAAAASASGSPGPGEGSAGGEKRSTAP
SAAASASASAAASSPAGGGAEALLEHCGVCRERLR
PEREPRLLPCLHSACSACLGPAAPAAANSSGDGGAA
GDGTVVDCPVCKQQCFSKDIVENYFMRDSGSKAATD
AQDANQCCTSCEDNAPATSYCVECSEPLCETCVEAH
QRVKYTKDHTVRSTGPAKSRDGERTVYCNVHKHEPL
VLFCESCDTLTCRDCQLNAHKDHQYQFLEDAVRNQR
KLLASLVKRLGDKHATLQKSTKEVRSSIRQVSDVQKRV
QVDVKMAILQIMKELNKRGRVLVNDQKVTEGQQERL
ERQHWMTMTKIQKHQEHILRFASWALESDNNTALLLSK
KLIYFQLHRALKMIVDPVEPHGEMKFQWDLNAWTKSA
EAFGKIVAERPNTSTGPAPMAPPRAPGPLSKQGSGS
SQPMEVQEGYGFGSGDDPYSSAEPHVSGVKRSRSG
EGEVSGLMRKVPRVSLERLDDLTDADSQPPVFKVFP
STTEDYNLIVIERGAAAAATGQPGTAPAGTPGAPPLAG
MAIVKEEETEAAGAPPTATEGPETKPVLMALAEPGA
EGPRLASPSGSTSSGLEVVAPGTSAPGGGPGTLDD
SATICRVCQKPGDLVMCNQCEFCFHLDCPLALQDVP
GEEWSCSLCHVLPDLKEEDGSLSLDGADSTGVVAKLS
PANQRKCEVLLALFCHEPCRPLHQLATDSTFSLDQP
GGTLDLTLIRARLQEKLSPPYSSPQEFQADVGRMFKQ
FNKLTEDKADVQSIIGLQRRFFETRMNEAFGDTKFS
AVLVEPPPMSLPGAGLSSQELSGGPGDGP

Host: Mouse

Reactivity: Human

Applications: IF, WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10155

Gene Symbol: TRIM28

Gene Alias: FLJ29029, KAP1, RNF96, TF1B, TIF1B

Gene Summary: The protein encoded by this gene mediates transcriptional control by interaction with the Kruppel-associated box repression domain found in many transcription factors. The protein localizes to the nucleus and is thought to associate with specific chromatin regions. The protein is a member of the tripartite motif family. This tripartite motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. [provided by RefSeq]