

Datasheet

ZNF256 monoclonal antibody (M02), clone 2H6

Catalog Number: H00010172-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant ZNF256.

Clone Name: 2H6

Immunogen: ZNF256 (NP_005764.2, 521 a.a. ~ 627
a.a) partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

CNECGKFFSQSSSLIRHRRSHTGERPYECSECWKSFS
NHSSLVKHRRVHTGERPYECSECGKSFSQSSNLTNH
QRIHSGERPYECSDCGKFFTFNSNLLKHQNVHKG

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10172

Gene Symbol: ZNF256

Gene Alias: BMZF-3, BMZF3