

Datasheet

AGR2 monoclonal antibody (M05), clone 1C9

Catalog Number: H00010551-M05

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant AGR2.

Clone Name: 1C9

Immunogen: AGR2 (AAH15503.1, 1 a.a. ~ 175 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MEKIPVSAFLLLVALSYTLARDTTVKPGAKKDTKDSRP
KLPQTL SRGWG DQLIWTQTYEEALYKSKTSNKPLMIH
HLDECPHSQALKKVFAENKEIQKLAEQFVLLNLVYETT
DKHLSPDGQYVPRIMFVDP SLTVRADITGRYSNRLYAY
EPADTALLLDNMKKALKLLKTEL

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re, WB-Ti, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10551

Gene Symbol: AGR2

Gene Alias: AG2, GOB-4, HAG-2, XAG-2