

Datasheet

RABAC1 monoclonal antibody (M01), clone 2A4

Catalog Number: H00010567-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant RABAC1.

Clone Name: 2A4

Immunogen: RABAC1 (NP_006414.1, 1 a.a. ~ 77 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MAAQKDQQKDAEAEGLSGTTLPLKIPSGAGREWLER
RRATIRPWSTFVDQQRFSRPRNLGELCQRLVRNVEYY
QSN

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10567

Gene Symbol: RABAC1

Gene Alias: PRA1, PRAF1, YIP3