

Datasheet

IFITM2 monoclonal antibody (M14), clone 1F2

Catalog Number: H00010581-M14

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant IFITM2.

Clone Name: 1F2

Immunogen: IFITM2 (NP_006426.1, 1 a.a. ~ 59 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MNHIVQTFSPVNSGQPPNYEMLKEEQEVAMLGAPHN
PAPPTSTVIHIRSETSVDPDHVWW

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Ce, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10581

Gene Symbol: IFITM2

Gene Alias: 1-8D