

Datasheet

HEXIM1 MaxPab mouse polyclonal antibody (B02)

muscle cells. This gene has no introns. [provided by RefSeq]

Catalog Number: H00010614-B02

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human HEXIM1 protein.

Immunogen: HEXIM1 (NP_006451.1, 1 a.a. ~ 359 a.a) full-length human protein.

Sequence:

```
MAEPFLSEYQHQPQTSNCTGAAAVQEELNPERPPGA  
EERVPEEDSRWQSRAFPQLGGRPGPEGEGLSLESQPP  
PLQTQACPESSCLREGEKGQNGDDSSAGGDFPPPAE  
VEPTPEAELLAQPCHDSEASKLGAPAAGGEEEWGQQ  
QRQLGKKKHRRRPSKKKRHWKPYKLTWEEKKKFDE  
KQSLRASRIRAEMFAKGQPVAPYNTTQFLMDDHDQEE  
PDLKTGLYSKRAAAKSDDTSDDDFMEEGGEEDGGSD  
GMGGDGSEFLQRDFSETYERYHTESLQNMSKQELIKE  
YLELEKCLSRMEDENNRLRLESKRLGGDDARVRELEL  
ELDRLRAENLQLLTENELHRQGERAPLSKFGD
```

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10614

Gene Symbol: HEXIM1

Gene Alias: CLP1, EDG1, FLJ13562, HIS1, MAQ1

Gene Summary: Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth