

Datasheet

ZNF263 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00010127-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human ZNF263 protein.

Immunogen: ZNF263 (NP_005732.2, 1 a.a. ~ 683 a.a) full-length human protein.

Sequence:

MASGPGSQEREGLLIVKLEEDCAWSQELPPDPGPSP
EASHLRFRFRFQEAAGPREALSRLQELCHGWL RPE
MRTKEQILELLVLEQFLTILPQEIQSRVQELHPESGEEA
VTLVEDMQRELGRRLRQQVTNHGRGTEVLLEEPLPLET
ARESPSFKLEPMETERSPGPRLQELLGPSPQRDPQAV
KERALSAPWLSLFPPEGNMEDKEMTGPQLPESLEDVA
MYISQEEWGHQDPSKRALS RDTVQESYENVDSLESHI
PSQEVPGTQVGQGGKLWDPSVQSCKEGLSPRG PAP
GEEKFENLEGVPSVCS ENIHPQVLLPDQARGEVPWSP
ELGRPHDRSQGDWAPPPEGGMEQALAGASSGRELG
RPKELQPKKLHLCPLCGKNFSNNSNLIRHQRIHAAERL
CMGVDCTEIFGGNPRFLSLHRAHLGEEAHKCLECGKC
FSQNTHLTRHQRTHTGEKPYQCNICGKCFSCNSNLHR
HQRTHTGEKPYKCPECGEIFAHSSNLLRHQRIHTGER
PYKCPEGKSFSSRSHLVIHERTHERERLYPFSECGE
AVSDSTPFLT NHGAHKA EKKLFECLTCGKSFRQGMHL
TRHQRTHTGEKPYKCTLCGENFSHRSNLIRHQRIHTG
EKPYTCHECGDSFSHSSNIRHLRHTHTGERPYKCSEC
GESFSRSSRLMSHQRTHTG

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to

avoid repeated freezing and thawing.

Entrez GeneID: 10127

Gene Symbol: ZNF263

Gene Alias: FPM315, ZKSCAN12