

Datasheet

CHST4 polyclonal antibody (A01)

Catalog Number: H00010164-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant CHST4.

Immunogen: CHST4 (NP_005760, 221 a.a. ~ 320 a.a)
partial recombinant protein with GST tag.

Sequence:

LMIDSRIVMGQHEQKLKKEDQPYVVMQVICQSQLEIYK
TIQSLPKALQERYLLVRYEDLARAPVAQTSRMYEFVGL
EFLPHLQTWVHNITRGKMGMDHAF

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10164

Gene Symbol: CHST4

Gene Alias: LSST